

















Exélia神话助听器

主要技术参数及功能对照表

技术参数		Exélia CIC/MC 标准版 /Petite 微型版	Exélia ITC/HS Petite 微型版	Exélia ITC/HS 标准版	Exélia FS	Exélia FS P	Exélia M	Exélia P	Exélia SP
									
最大声输出	2cc耦合腔	113	113	113	115	122	125	130	131
	仿真耳	123	123	123	125	130	131	138	139
满档增益	2cc耦合腔	35	48	48	50	55	53	63	70
	仿真耳	45	59	59	60	63	58	68	75
频响范围	2cc耦合腔	<100 - 7000	<100 - 6600	<100 - 6600	<100 - 6400	<100 - 5500	<100 - 6800	<100 - 6500	<100 - 6500
	仿真耳	150 - 7000	150 - 6800	150 - 6800	100 - 6600	100 - 5600	100 - 7000	100 - 7000	100 - 6800
电池	型号	10A	312	312GenX	13GenX	13GenX	13GenX	13GenX	13GenX
	耗电量	1.2	1.2	1.2	1.2	1.2	1.3	1.3	1.3
验配范围									
多通道自适应言语聚焦功能			✓	✓	✓	✓	✓	✓	✓
真耳声音重塑							✓	✓	✓
一键式双耳同步遥控				✓	✓	✓	✓	✓	✓
言语声源方位遥控功能				✓	✓	✓	✓	✓	✓
可视遥控器myPilot				✓	✓	✓	✓	✓	✓
蓝牙/音频转换器iCom				✓	✓	✓	✓	✓	✓
FM			✓	✓	✓	✓	✓	✓	✓

注：①Petite微型版定制机指为减小体积而不带无线互通和遥控功能的定制机。

②GenX为高性能助听器电池，为确保助听器性能，所有带无线互通功能和遥控功能的助听器均需使用GenX专用电池。